

UCS[®] DEBRIDEMENT

Sterile system for initial and maintenance wound debridement, cleansing and hydration of the peri-lesional area as well as the whole limb

PRODUCT FEATURES:

Atraumatic removal of:

- Fibrin.
- Necrotic tissue.
- Slough.
- Hyperkeratotic edge.
- Bacterial bio-film.

The UCS[®] Debridement gauze captures and retains the debris:

- Cleansing and hydration of the peri-lesional area as well as the whole limb.
- No bleeding and damage to the granulation tissue.
- Reduces pain, increasing patient compliance.
- Shortens the time required for the procedure.
- The gentle, yet effective debridement allows the use of the product at every dressing change.
- Recommended for the treatment of venous, arterial, mixed, diabetic, decubitus (I-IV stage) abscesses and removal of residues of previous dressings.
- **UCS[®] Debridement** is used in major hospitals and university centers in Europe.



AREA OF USE:

- Chronic wounds.
- Acute wounds.
- 1° and 2° degree burns.
- Fistulas and abscesses.
- Peristomal skin.
- Ports of entry of catheters, PEG/PEJ.
- Removal of bandages encrusted to the wound.

INGREDIENTS:

Aqua

Aloe barbadensis leaf extract

Allantoin

Poloxamer

TECHINCAL FEATURES:

Viscose to provide:

- High absorbency of skin impurities.
- High resistance to stretch and torsion to avoid tearing of the wipe.
- High softness to avoid irritations.
- High solution content in each gauze. (12 g)
- Gauze size (open) 20 cm x 20 cm
- Sterile.

CLASSIFICATION:

Sterile Class IIb Medical Device.

CERTIFICATIONS AND ENVIRONMENTAL CONTROL:

Welcare is certified:

UNI CEI EN ISO 13485 • UNI EN ISO 9001 • UNI EN ISO 14001

Productive process:

Manufactured in clean room



Bibliography:

- ♦ Efficacy of anti bacterial wipes [The Institute for skin research 11.05.2008]
- ♦ 5 Day test for MRSA on pig skin [ITS Testing Service Ltd 28.11.2008]
- ♦ In vivo and in vitro studies for the anti bacterial wipes efficacy [The institute for skin research 11.11.2008]
- ♦ Cytotoxicity direct contact [Eurofins Pharma Services - 24.06.2009]
- ♦ A new Ulcer Cleansing System in the management of "no healing" ulcers in elderly. [Palumbo, Serantoni, Abbritti, Falla, Di Salvo, 17-18/04/2013]
- ♦ Faulkner DM, Sutton ST, Hesford JD, Faulkner BC, Major DA, Hellewell TB, et al., A new stable pluronic F68 gel carrier for antibiotics in contaminated wound treatment. Am J Emerg Med. 1997; 15(1):
- ♦ Bryant CA, Rodeheaver GT, Reem EM, Nichter LS, Kenney JG, Edlich RF. Search for a nontoxic surgical scrub solution for periorbital lacerations. An Emerg Med. May 1984; 13(5):317 -321.
- ♦ Rodeheaver GT, Kurtz L, Kircher BJ, Edlich RF, Pluronic F-68: a promising new skin wound cleanser. Ann Emerg Med. 1988; 9(11): 572 -6.
- ♦ Curry DJ, Wright DA, Lee RC, Kang UJ, Frim DM. Surfactant poloxamer 188 - related decreases in inflammation and tissue damage after experimental brain injury in rats. J Neurosurg. Aug 2004;101(1 Suppl):91 -96.

- ♦ Birchenough, SA, MD; Peirce, SM, PhD; Rodeheaver, G, PhD; Morgan, RF, MD; Katz, AJ, MD. Topical Poloxamer-188 Improves Blood Flow Following Thermal Injury. Abstracted presented at the Southeastern Society of Plastic and Reconstructive Surgery Conference [SESPRS] June 2007.
- ♦ Mustafi D, Smith CM, Makinen MW, Lee RC. Multi -block poloxamer surfactants suppress aggregation of denatured proteins. Biochim Biophys Acta. 2008;1780(1):7 -15.
- ♦ Tursilli R, Piel G, Delattre L, Scalia S. Solid lipid microparticles containing the sunscreens agent, octyl-dimethylaminoazobenzoate: effect of the vehicle. Eur J. Pharm Biopharm. Jun 2007;66(3):483 -487.
- ♦ A new dimension in wound cleansing [R. Cassino - 19th EWMA Conference 20-22 May 2009 - Helsinki Finland]
- ♦ La deterzione come sistema di rimozione dello stafilococco aureo meticillinoresistente (MRSA) dalle lesioni cutanee non infette [R. Cassino - VIII Congresso Nazionale AIUC 23 -26 Settembre 2009 - Firenze]
- ♦ Efficacy and tolerability evaluation of UCS Wound Cleansing System. Clinical study [Pori City Hospital - Finland 2009]
- ♦ Efficacy and tolerability evaluation of UCSol Wound Cleansing Solution vs. other ringier solutions and saline solution [Homecare Service Florence - Italy 2009]
- ♦ A new ulcer cleansing system in the management of "non healing" ulcers in elderly [F.Palumbo, S.Serantoni, F.Abbritti, G.Falla, M.M.Di Salvo. Tissue Viability Society, Kittering, U.K. 17th -18th April 2013.
- ♦ Evaluation of a new mechanical debridement system [Serantoni S., Palumbo F.P., Abbritti F. 2012]

- ♦ Effective debridement can be achieved in a busy clinic environment [Audrey Gillies, medi uk, 2015]
- ♦ Cost efficacy of using Juxtra CURESTM and UCSTM debridement cloths [Sue Elvin, nurse consultant, district nursing, Camden IPC, 2015]
- ♦ How one trust saved money by changing its leg cleansing methods in a switch to debridement cloths [Maria A Hughes, Tissue Viability Lead Specialist Nurse, Wirral Community NHS Trust, Birkenhead, Wirral, 2015]
- ♦ Changing clinical practice in leg ulcer management, implementation of a new leg ulcer cleansing system [Maria A Hughes, Tissue Viability Lead Specialist Nurse, Wirral Community NHS Trust, Birkenhead, Wirral, 2015]
- ♦ How wound cleansing and debriding aids management and healing [Annette Downe, clinical nurse specialist, St. Thomas Hospital, London]
- ♦ Ma.De.U.L. Study [Maintenance Debridement Ulcer and Leg Study] Giovanni Mosti, Stefano Gasperini Angiology Department, M.D. Barbantini Clinic - Lucca (Italy) Medical Advisor - Pisa (Italy)
- ♦ Traumatic foot injury with skin graft failure back to basic wound care fundamentals [Rene Lessing, Rn, Dipl Nursing, Community Nursing, Psychiatry, Midwifery, Vkovs, Cert. Woundcare Uf, Cert. Woundcare University Of Hertfordshire - United Kingdom.]

CODE:

DESCRIPTION:

SALES UNITS:

UCS® Debridement

95150

1 pouch: 1 pre-moistened wipe. 1 box: 30 pouches

1 carton: 10 boxes

95150F

1 pouch: 1 pre-moistened wipe. 1 box: 10 pouches

1 carton: 30 boxes

Welcare Industries S.p.A.

Operational site: via dei Falegnami, 7 - 05018 Orvieto (TR)

Commercial office: Regus Building, piazzale Biccamano, 8 - 20121 Milano

Commercial office: EUR - Spaces, viale Luca Gaurico, 91/93 - 00143 Roma

Registered office: via San Giovanni sul Muro, 18 - 20121 Milano

www.welcareindustries.com

Revision 00 2022/03

WELCARE
Protecting People - Caring for Nature